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
Power Output

Users of old style vertical arc lamp housings are throwing away as much as 90% of the lamps output, due to poor collection efficiency. These old style vertical lamp housings have a collection lens in front of the arc lamp and sometimes, but not always, a back reflector behind them. The problem with this old design is that only the light that actually strikes these optical elements is delivered outside of the lamp housing. All other photons emitted by the lamp are wasted, simply heating the inside of the lamp housing. Conversely the unique OBB PowerArc™ lamp housing has an enveloping ellipsoidal reflector that collects virtually all of the light emitted by the lamp arc, delivering those photons to a secondary focal point outside of the lamp housing, and it does so without any lenses.

For the same bulb, the OBB PowerArc™ delivers 5 to 6 times more light, to the secondary focus! That means that an OBB PowerArc™ lamp housing with a 75 watt xenon lamp provides the equivalent optical power of a 450 watt xenon lamp in an old style vertical lamp housing. And it does so with greater power density due to a smaller focal spot, and at a small fraction of the cost of a big old 450 watt illuminator.

Lamp Housing

At the heart of the OBB PowerArc™ lamp housing is a proprietary on-axis ellipsoidal reflector. Our reflectors collect up to 70% of the radiant energy from the arc lamp, versus only 12% for typical condenser systems in vertical lamp housings. The ellipse literally wraps around the arc lamp, collecting 5 to 6 times more output power than a conventional system.


The arc source is located at one focal point of the ellipse, and the radiation is reflected by the ellipse to the secondary focus which is actually outside of the lamp housing. Since the light is brought to a focus by reflection rather than refraction (through a lens), there are less losses from absorption or lens-surface back-reflection. This design is so efficient that an OBB PowerArc™ lamp housing can deliver up to 11 times more optical power into a given smaller area than a conventional lamp housing. This is critical when illuminating light guides, monochromator slits, pinholes or other small areas.

What this means is simply that you get the same output with an OBB 75 W system as with a conventional 450 W system. You obviously will save money and space.

While conventional lamp housings resemble chimneys emitting ozone and requiring cumbersome venting, the OBB PowerArc™ has a sealed lamp housing that requires no ozone venting.

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